

## Safety Considerations

### Recommendations

Dräger recommends the use of an Uninterruptible Power Manager in conjunction with the Infinity CentralStation at all times.

### General Precautions

**WARNING! Because of the danger of electric shock, never remove the cover of any device while it is in operation or connected to a power outlet.**

The user must check that the equipment functions safely and is in proper working condition before use. In the interest of safety, regular equipment inspection and maintenance is required. Once a year, check all cables, devices, and accessories for damage, ground resistance, chassis and patient leakage currents, and all alarm functions. Also, ensure that all safety labels are legible. Maintain a record of safety inspections. Safety checks, device verification, calibration, and maintenance should be carried out and verified by Dräger-authorized personnel.

#### **WARNING!**

- **Operation of this device should be in accordance with part 15.242 sections A, D, E, F, G and H of the FCC rules.**
- **Take all precautions against electric shock due to exposed electrical connections, spilled fluids, defibrillation, etc.**

**CAUTION!** *Conductive parts of electrodes and associated connectors for applied parts, including the neutral electrode, should not contact other conductive parts including earth.*

**NOTE:** This equipment has been tested and found to comply with the limits for a Class B digital device, pursuant to Part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference in a residential installation. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instructions, may cause harmful interference to radio communications. However, there is no guarantee that interference will not occur in a particular installation. If this equipment does cause harmful interference to radio or television reception, which can be determined by turning the equipment off and on, the user is encouraged to try to correct the interference by one or more of the following measures:

- Reorient or relocate the receiving antenna,
- Increase the separation between the equipment and receiver,
- Connect the equipment into an outlet on a different circuit from that to which the receiver is connected,
- Consult the dealer or an experienced radio/TV technician for help.

*Changes or modifications to this equipment not expressly approved by the party responsible for compliance could void the user's authority to operate the equipment.*